




Product specification

Mailing envelope, PE, (LxW): 310 x 180 mm



Product information

Product number	95.1410
----------------	---------

Product descriptions

Name	Mailing envelope, PE, (LxW): 310 x 180 mm
Product number	95.1410
Product description	Mailing envelope, material: PE, blue, (LxW): 310 x 180 mm, 50 piece(s)/case

Material & colours

Color of product	blue
------------------	------

Size

Width of product	180 mm
Length of product	310 mm

Product characteristics

Suction insert	no
----------------	----

Purity & certification

Satisfies the requirement	P 650
Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / inner box	50
Piece(s) / case	1000
Piece(s) / pallet	64000
Height of inner box	20 mm
Width of inner box	190 mm
Depth of inner box	320 mm
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	11.5 kg
EAN inner box	4038917086530
EAN case	4038917005081

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT