



Product specification

Rack, styrofoam, format: 20 x 5, suitable for S-Monovette® Ø 11 mm



Product information

Product number	95.1439
----------------	---------

Product descriptions

Name	Rack, styrofoam, format: 20 x 5, suitable for S-Monovette® Ø 11 mm
Product number	95.1439
Product description	Rack, material: styrofoam, white, format: 20 x 5, (LxWxH): 338 x 90 x 50 mm, suitable for S-Monovette® Ø 11 mm, 1 piece(s)/blister

Material & colours

Color of product	white
------------------	-------

Size

Width of product	90 mm
Length of product	338 mm
Height of product	50 mm

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	40
Piece(s) / pallet	640
Height of inner box	53 mm
Width of inner box	85 mm
Depth of inner box	310 mm
Height of case	345 mm
Width of case	405 mm
Depth of case	455 mm
Volume of case	0.0636 cbm
Weight of case	1.7 kg
EAN inner box	4038917086646
EAN case	4038917005173

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT