



Product specification

Transport case B 17, with lateral rim for block rack



Product information

Product number	95.1717
----------------	---------

Product descriptions

Name	Transport case B 17, with lateral rim for block rack
------	------------------------------------------------------

Product number	95.1717
----------------	---------

Product description	Transport case B 17, with lateral rim for block rack, (LxWxH): 395 x 395 x 435 mm
---------------------	-----------------------------------------------------------------------------------

Size

Width of product	395 mm
------------------	--------

Length of product	395 mm
-------------------	--------

Height of product	435 mm
-------------------	--------

Product characteristics

Design	with lateral rim for block rack
--------	---------------------------------

Purity & certification

Satisfies the requirement	P 650, ADR
---------------------------	------------

Product category	no medical device IVD
------------------	-------------------------

Batched	no
---------	----

Piece(s) / pallet	18
Height of inner box	123 mm
Width of inner box	155 mm
Depth of inner box	352 mm
Height of case	460 mm
Width of case	403 mm
Depth of case	403 mm
Volume of case	0.0747 cbm
Weight of case	4.5 kg
EAN inner box	4038917087544
EAN case	4038917065610

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT