



## Product specification

### Post transport packaging, 198 x 107 x 38 mm, for diagnostic specimens



## Product information

Product number	95.900
----------------	--------

## Product descriptions

Name	Post transport packaging, 198 x 107 x 38 mm, for diagnostic specimens
------	---

Product number	95.900
----------------	--------

Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, flat design, (LxWxH): 198 x 107 x 38 mm, for diagnostic specimens, holds 1 to 4 transport containers/bottles up to a diameter of 30 mm, Locking mechanism, 50 piece(s)/case
---------------------	---

## Material & colours

Color of product	brown
------------------	-------

## Size

Width of product	107 mm
------------------	--------

Length of product	198 mm
-------------------	--------

Height of product	38 mm
-------------------	-------

Design	flat design
Label/Print	With print
Color of print/label	black
Example of use	External packaging according to P 650/ ADR

## Purity & certification

Satisfies the requirement	P 650, ADR, IATA
Product category	not regulated
CE certified	no CE mark
Batched	yes

## Packaging

Piece(s) / inner box	50
Piece(s) / case	50
Piece(s) / pallet	3000
Height of case	115 mm
Width of case	366 mm
Depth of case	584 mm
Volume of case	0.0246 cbm
Weight of case	3.44 kg
EAN inner box	4038917084550
EAN case	4038917064576

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT