



## Product specification

**Post transport packaging, 107 x 198 x 50 mm, for diagnostic specimens**



## Product information

Product number	95.901
----------------	--------

## Product descriptions

Name	Post transport packaging, 107 x 198 x 50 mm, for diagnostic specimens
Product number	95.901
Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, high design, (LxWxH): 107 x 198 x 50 mm, for diagnostic specimens, holds 1 to 4 transport containers/bottles up to a diameter of 30 mm, Locking mechanism, 50 piece(s)/case

## Material & colours

Product material	Cardboard
Color of product	brown

## Size

Width of product	198 mm
Length of product	107 mm
Height of product	50 mm

## Product characteristics

Design	high design
Label/Print	With print
Color of print/label	black
Example of use	External packaging according to P 650/ ADR

## Purity & certification

Satisfies the requirement	P 650, ADR, IATA
Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	50
Piece(s) / inner box	50
Piece(s) / case	50
Piece(s) / pallet	2100
Height of case	115 mm
Width of case	414 mm
Depth of case	656 mm
Volume of case	0.0312 cbm
Weight of case	4.04 kg
EAN inner box	4038917084567
EAN case	4038917064583

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT