



Product specification

Post transport packaging, 220 x 170 x 40 mm, for diagnostic specimens



Product information

Product number	95.902
----------------	--------

Product descriptions

Name	Post transport packaging, 220 x 170 x 40 mm, for diagnostic specimens
Product number	95.902
Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, flat design, (LxWxH): 220 x 170 x 40 mm, for diagnostic specimens, holds 1 to 4 transport containers/bottles up to a diameter of 30 mm, enclosed foam foil pouch, 50 piece(s)/case

Material & colours

Color of product	brown
------------------	-------

Size

Width of product	170 mm
Length of product	220 mm
Height of product	40 mm

Design	flat design
Label/Print	With print
Color of print/label	black
Example of use	External packaging according to P 650/ ADR

Purity & certification

Satisfies the requirement	P 650, ADR, IATA
Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	50
Piece(s) / pallet	2400
Height of case	145 mm
Width of case	401 mm
Depth of case	488 mm
Volume of case	0.0284 cbm
Weight of case	3.86 kg
EAN inner box	4038917084574
EAN case	4038917069854

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT