



Product specification

Felt marker, blue, waterproof



Product information

Product number	95.953
----------------	--------

Product descriptions

Name	Felt marker, blue, waterproof
------	-------------------------------

Product number	95.953
----------------	--------

Product description	Felt marker, blue, waterproof, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case
---------------------	---

Material & colours

Color of product	black
------------------	-------

Size

Diameter	10 mm
----------	-------

Length of product	142 mm
-------------------	--------

Product characteristics

Color of print/label	white
----------------------	-------

Example of use	for labelling plastic tubes and reaction tubes
----------------	--

Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
--------------	------------

Batched	yes
---------	-----

Packaging

Piece(s) / inner box	10
Piece(s) / case	1000
Piece(s) / pallet	48000
Height of inner box	28 mm
Width of inner box	70 mm
Depth of inner box	150 mm
Height of case	292 mm
Width of case	300 mm
Depth of case	392 mm
Volume of case	0.0343 cbm
Weight of case	10 kg
EAN inner box	4038917084628
EAN case	4038917003391

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT