



# Product specification

## Felt marker, black, waterproof



### Product information

Product number	95.954
----------------	--------

### Product descriptions

Name	Felt marker, black, waterproof
Product number	95.954
Product description	Felt marker, black, waterproof, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case

### Material & colours

Product material	Polypropylene (PP)
Color of product	black

### Size

Diameter	10 mm
Length of product	142 mm

### Product characteristics

Color of print/label	white
Example of use	for labelling plastic tubes and reaction tubes

### Purity & certification

Product category	no medical device   IVD
------------------	-------------------------

Batched

yes

## Packaging

Minimum order quantity	10
Piece(s) / inner box	10
Piece(s) / case	1000
Piece(s) / pallet	48000
Height of case	392 mm
Width of case	291 mm
Depth of case	300 mm
Volume of case	0.0342 cbm
Weight of case	10 kg
EAN inner box	4038917084635
EAN case	4038917003407

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT