



Product specification

IsoFreeze® PCR Rack, PP, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, PCR strips and PC



Product information

Product number	95.984
----------------	--------

Product descriptions

Name	IsoFreeze® PCR Rack, PP, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, PCR strips and PCR plates
Product number	95.984
Product description	IsoFreeze® PCR Rack, material: PP, multi-coloured, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, PCR strips and PCR plates, Colour change from purple to pink at 7 °C, 2 piece(s)/case

Material & colours

Color of product	multi-coloured
Color of frame	white
Color of cap	transparent

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	10
Piece(s) / pallet	1280
Height of inner box	181 mm
Width of inner box	36 mm
Depth of inner box	181 mm
Height of case	175 mm
Width of case	197 mm
Depth of case	584 mm
Volume of case	0.0201 cbm
Weight of case	2.95 kg
EAN inner box	4038917373425
EAN case	4038917373418

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT