




## Product specification

### Base station, PC, suitable for work tray 95.987.002 & 0.2 ml PCR tubes



## Product information

Product number	95.988
----------------	--------

## Product descriptions

Name	Base station, PC, suitable for work tray 95.987.002 & 0.2 ml PCR tubes
Product number	95.988
Product description	Base station, material: PC, transparent, (LxWxH): 128 x 85.5 x 30.4 mm, suitable for work tray 95.987.002 & 0.2 ml PCR tubes, 5 piece(s)/case

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Width of product	85.5 mm
Length of product	128 mm
Height of product	30.4 mm

## Purity & certification

Product category	not regulated
CE certified	no CE mark

Batched	no
---------	----

## Packaging

Piece(s) / case	75
Piece(s) / pallet	3600
Height of case	145 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0345 cbm
Weight of case	2.83 kg
EAN inner box	4038917084864
EAN case	4038917003636

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT