

Product specification

Pressure infusion cuff K



Product information

Order number	97.8100.000
--------------	-------------

Product descriptions

Name	Pressure infusion cuff K
Product number	97.8100.000
Product description	Pressure infusion cuff K, up to 1000 ml, variable fit thanks to velcro fastener

Product characteristics

Type of cap	Velcro fastener
-------------	-----------------

Size

Width of product	16 cm
Length of product	28 cm
Height of product	8 mm
Max. volume	1000 ml

Material & colours

Color of product	red
------------------	-----

Purity & certification

Product category	Medical device
CE certified	CE
Latex-free	no
Batched	no

Packaging

Volume of case	0.0039 cbm
Weight of case	0.52 kg
EAN inner box	4028339310122

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

