

Product specification

Pressure infusion cuff Standard



Product information

Order number	97.8100.100
--------------	-------------

Product descriptions

Name	Pressure infusion cuff Standard
------	---------------------------------

Product number	97.8100.100
----------------	-------------

Product description	Pressure infusion cuff Standard, up to 500 ml, fixed insertion
---------------------	--

Product characteristics

Type of cap	fixed insertion
-------------	-----------------

Size

Width of product	16 cm
------------------	-------

Length of product	28 cm
-------------------	-------

Height of product	8 cm
-------------------	------

Max. volume	500 ml
-------------	--------

Material & colours

Color of product	red
------------------	-----

Purity & certification

Accordance with MDD or IVDD	Medical device
-----------------------------	----------------

CE certified	CE
--------------	----

Latex-free	no
------------	----

Lot or serial number	no
----------------------	----

Packaging

Volume of case	0.0039 cbm
----------------	------------

Weight of case	0.39 kg
----------------	---------

EAN inner box	4028339310825
---------------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

