




Product specification

Double block rack D13, Ø opening: 13 mm, 5 x 20, grey



Product information

Product number	93.853.300
----------------	------------

Product descriptions

Name	Double block rack D13, Ø opening: 13 mm, 5 x 20, grey
Product number	93.853.300
Product description	Double block rack D13, for 100 collection tubes, Ø opening: 13 mm, format: 5 x 20, grey, material: PP

Material & colours

Product material	Polypropylene (PP)
Color of product	grey

Size

Width of product	96 mm
Length of product	360 mm
Height of product	45 mm
Diameter of opening	13 mm

Product characteristics

Number of wells	100
Storage capacity	100
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / pallet	720
Height of case	120 mm
Width of case	368 mm
Depth of case	398 mm
Volume of case	0.0176 cbm
Weight of case	3.27 kg
EAN inner box	4038917122634
EAN case	4038917074087

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT