




Product specification

Rack, styrofoam, format: 10 x 5, suitable for 15 ml tubes



Product information

Product number	95.64.251
----------------	-----------

Product descriptions

Name	Rack, styrofoam, format: 10 x 5, suitable for 15 ml tubes
Product number	95.64.251
Product description	Rack, material: styrofoam, white, format: 10 x 5, (LxWxH): 205 x 127 x 46 mm, for 50 collection tubes, suitable for 15 ml tubes, 1 piece(s)/blister

Material & colours

Product material	Styrofoam
Color of product	white
Closure material	Styrofoam
Color of cap	white

Size

Width of product	127 mm
Length of product	205 mm
Height of product	46 mm

Product characteristics

Storage capacity	50
------------------	----

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	21
Piece(s) / case	21
Piece(s) / pallet	882
Height of case	235 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0345 cbm
Weight of case	0.58 kg
EAN inner box	4038917095594
EAN case	4038917007771

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT