




Product specification

Rack, styrofoam, format: 20 x 5, suitable for tubes Ø 15 and 16 mm



Product information

Product number	95.64.250
----------------	-----------

Product descriptions

Name	Rack, styrofoam, format: 20 x 5, suitable for tubes Ø 15 and 16 mm
Product number	95.64.250
Product description	Rack, material: styrofoam, white, format: 20 x 5, (LxWxH): 416 x 107 x 50 mm, for 100 collection tubes, suitable for tubes Ø 15 and 16 mm, 1 piece(s)/blister

Material & colours

Product material	Styrofoam
Color of product	white
Closure material	Styrofoam
Color of cap	white

Size

Width of product	107 mm
Length of product	416 mm
Height of product	50 mm

Product characteristics

Storage capacity	100
------------------	-----

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / case	28
Piece(s) / pallet	448
Height of case	380 mm
Width of case	428 mm
Depth of case	498 mm
Volume of case	0.081 cbm
Weight of case	1.65 kg
EAN inner box	4038917095587
EAN case	4038917007764

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT