



Product specification

IsoFreeze® MCT Rack, PP, format: 6 x 4, suitable for reaction tubes Ø 10.8 mm (1.5 ml & 2 ml)



Product information

Product number	95.983
----------------	--------

Product descriptions

Name	IsoFreeze® MCT Rack, PP, format: 6 x 4, suitable for reaction tubes Ø 10.8 mm (1.5 ml & 2 ml)
Product number	95.983
Product description	IsoFreeze® MCT Rack, material: PP, multi-coloured, format: 6 x 4, suitable for reaction tubes Ø 10.8 mm (1.5 ml & 2 ml), 1 piece(s)/case

Material & colours

Color of product	multi-coloured
Color of frame	white
Color of cap	transparent

Product characteristics

Label/Print	With print
-------------	------------

Purity & certification

Product category	not regulated
CE certified	no CE mark

Batched	no
---------	----

Packaging

Piece(s) / pallet	240
Height of inner box	95 mm
Width of inner box	115 mm
Depth of inner box	190 mm
Height of case	93 mm
Width of case	197 mm
Depth of case	584 mm
Volume of case	0.0107 cbm
Weight of case	2.5 kg
EAN inner box	4038917394895
EAN case	4038917394888

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT