




Product specification

Rack, PC, format: 10 x 2, suitable for tubes, all S-Monovette® diameters



Product information

Product number	93.1097.100
----------------	-------------

Product descriptions

Name	Rack, PC, format: 10 x 2, suitable for tubes, all S-Monovette® diameters
Product number	93.1097.100
Product description	Rack, material: PC, transparent, format: 10 x 2, (LxWxH): 257 x 74 x 55 mm, for 20 collection tubes, suitable for tubes, all S-Monovette® diameters, 1 piece(s)/case

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	74 mm
Length of product	257 mm
Height of product	55 mm

Product characteristics

Storage capacity	20
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	30
Piece(s) / pallet	900
Height of case	340 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0539 cbm
Weight of case	7.27 kg
EAN inner box	4038917125512
EAN case	4038917020121

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT