



Product specification

Rack, PP, format: 6 x 3, suitable for 50 ml centrifuge tubes



Product information

Product number	93.1707
----------------	---------

Product descriptions

Name	Rack, PP, format: 6 x 3, suitable for 50 ml centrifuge tubes
Product number	93.1707
Product description	Rack, material: PP, blue, format: 6 x 3, (LxWxH): 295 x 114 x 65 mm, for 18 collection tubes, suitable for 50 ml centrifuge tubes, 1 piece(s)/bag

Material & colours

Color of product	blue
------------------	------

Size

Width of product	114 mm
Length of product	295 mm
Height of product	65 mm

Product characteristics

Storage capacity	18
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Packaging

Piece(s) / case	10
Piece(s) / pallet	540
Height of inner box	105 mm
Width of inner box	130 mm
Depth of inner box	200 mm
Height of case	75 mm
Width of case	120 mm
Depth of case	300 mm
Volume of case	0.0027 cbm
Weight of case	2.81 kg
EAN inner box	4038917347501
EAN case	4038917347495

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT