



Product specification

Rack, PP, format: 6 x 3, suitable for 50 ml centrifuge tubes



Product information

Product number	93.1707
----------------	---------

Product descriptions

Name	Rack, PP, format: 6 x 3, suitable for 50 ml centrifuge tubes
Product number	93.1707
Product description	Rack, material: PP, blue, format: 6 x 3, (LxWxH): 295 x 114 x 65 mm, for 18 collection tubes, suitable for 50 ml centrifuge tubes, 1 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	blue

Size

Width of product	114 mm
Length of product	295 mm
Height of product	65 mm

Product characteristics

Storage capacity	18
Autoclavable (max. °C)	121 °C

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / case	10
Piece(s) / pallet	540
Height of case	300 mm
Width of case	75 mm
Depth of case	125 mm
Volume of case	0.0028 cbm
Weight of case	2.81 kg
EAN inner box	4038917347501
EAN case	4038917347495

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT