




Product specification

Protective cover, PC, brown, 256 x 62 x 72 mm, suitable for stand no 93.1097.100



Product information

Product number	93.1102.001
----------------	-------------

Product descriptions

Name	Protective cover, PC, brown, 256 x 62 x 72 mm, suitable for stand no 93.1097.100
Product number	93.1102.001
Product description	Protective cover, material: PC, brown, (LxWxH): 256 x 62 x 72 mm, suitable for stand no 93.1097.100, 1 piece(s)/bag

Material & colours

Product material	Polycarbonate (PC)
Color of product	brown

Size

Width of product	62 mm
Length of product	256 mm
Height of product	72 mm

Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / case	30
Piece(s) / pallet	1080
Height of case	300 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0475 cbm
Weight of case	4.61 kg
EAN inner box	4038917125529
EAN case	4038917020138

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT