



Product specification

Rack, PC, format: 6 x 3, suitable for tubes, square cuvettes, all S-Monovette® diameters



Product information

| | |
|----------------|------------|
| Product number | 93.850.001 |
|----------------|------------|

Product descriptions

| | |
|---------------------|--|
| Name | Rack, PC, format: 6 x 3, suitable for tubes, square cuvettes, all S-Monovette® diameters |
| Product number | 93.850.001 |
| Product description | Rack, material: PC, transparent, format: 6 x 3, (LxWxH): 137 x 70 x 40 mm, for 18 collection tubes, suitable for tubes, square cuvettes, all S-Monovette® diameters, 1 piece(s)/case |

Material & colours

| | |
|------------------|-------------|
| Color of product | transparent |
|------------------|-------------|

Size

| | |
|-------------------|--------|
| Width of product | 70 mm |
| Length of product | 137 mm |
| Height of product | 40 mm |

Product characteristics

| | |
|------------------|----|
| Storage capacity | 18 |
|------------------|----|

Autoclavable (max. °C)

121 °C

Purity & certification

| | |
|------------------|---------------|
| Product category | not regulated |
| CE certified | no CE mark |
| Batched | no |

Packaging

| | |
|---------------------|---------------|
| Piece(s) / case | 36 |
| Piece(s) / pallet | 2376 |
| Height of inner box | 138 mm |
| Width of inner box | 42 mm |
| Depth of inner box | 72 mm |
| Height of case | 155 mm |
| Width of case | 268 mm |
| Depth of case | 438 mm |
| Volume of case | 0.0182 cbm |
| Weight of case | 4.55 kg |
| EAN inner box | 4038917121750 |
| EAN case | 4038917018005 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT