



## Product specification

**Urine container, 3 l, with inspection strip and enclosed label with instructions for use, brown, wit**



## Product information

Product number	77.575.001
----------------	------------

## Product descriptions

Name	Urine container, 3 l, with inspection strip and enclosed label with instructions for use, brown, with light protection, graduated
Product number	77.575.001
Product description	Urine container, Urine collection, max. working volume: 3 l, Ø opening: 83.5 mm, brown, with light protection, graduated, material: PE, with paper label, english, cap: yellow, with inspection strip and enclosed label with instructions for use, 30 piece(s)/case

## Material & colours

Color of product	brown
Closure material	Polypropylene (PP)
Color of cap	yellow

## Size

Length excluding cap	243 mm
Length including cap	246 mm

Width of product	130 mm
Length of product	125 mm
Diameter of opening	83.5 mm
Max. volume	3 l
Sample volume	3000 ml

## Product characteristics

Design	with inspection strip and enclosed label with instructions for use
Label/Print	With paper label
Graduation	yes
Light protection	yes
Viewing window	yes
Closure type	Screw cap

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes
Safety data sheet	no

## Packaging

Piece(s) / inner box	30
Piece(s) / case	30
Piece(s) / pallet	270
Height of inner box	515 mm
Width of inner box	403 mm
Depth of inner box	643 mm
Height of case	515 mm
Width of case	403 mm
Depth of case	643 mm
Volume of case	0.1335 cbm
Weight of case	4.21 kg
EAN inner box	4038917266222
EAN case	4038917266215

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT