




## Product specification

**Pipette tip, 1,250 µl, transparent, Biosphere® plus, 96 piece(s)/SingleRefill**



## Product information

Product number	70.3060.305
----------------	-------------

## Product descriptions

Name	Pipette tip, 1,250 µl, transparent, Biosphere® plus, 96 piece(s)/SingleRefill
Product number	70.3060.305
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/SingleRefill

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	9 mm
Length of product	90 mm
Max. volume	1250 µl

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	69120
Height of inner box	148 mm
Width of inner box	105 mm
Depth of inner box	495 mm
Height of case	237 mm
Width of case	325 mm
Depth of case	512 mm
Volume of case	0.0394 cbm
Weight of case	4.3795 kg
EAN inner box	4038917432474
EAN case	4038917432467

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT