




## Product specification

**Filter tip, XL, 1,000 µl, transparent, Biosphere® plus, 96 piece(s)/box**



## Product information

Product number	70.3060.255
----------------	-------------

## Product descriptions

Name	Filter tip, XL, 1,000 µl, transparent, Biosphere® plus, 96 piece(s)/box
Product number	70.3060.255
Product description	Filter tip, working volume: 1,000 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, XL, 96 piece(s)/box

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	9 mm
Length of product	90 mm
Max. volume	1000 µl

## Product characteristics

Filter	yes
Graduation	yes

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	61440
Height of inner box	148 mm
Width of inner box	111 mm
Depth of inner box	530 mm
Height of case	245 mm
Width of case	325 mm
Depth of case	547 mm
Volume of case	0.0436 cbm
Weight of case	5.4106 kg
EAN inner box	4038917432269
EAN case	4038917432252

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT