




Product specification

Pipette tip, 1,250 µl, transparent, PCR Performance

Tested, Low retention, 96 piece(s)/box



Product information

Product number	70.3060.210
----------------	-------------

Product descriptions

Name	Pipette tip, 1,250 µl, transparent, PCR Performance Tested, Low retention, 96 piece(s)/box
Product number	70.3060.210
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	9 mm
Length of product	90 mm
Max. volume	1250 µl

Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	1920
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	61440
Height of inner box	148 mm
Width of inner box	111 mm
Depth of inner box	530 mm
Height of case	245 mm
Width of case	325 mm
Depth of case	547 mm
Volume of case	0.0436 cbm
Weight of case	5.2301 kg
EAN inner box	4038917432351
EAN case	4038917432344

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT