




## Product specification

**Pipette tip, 1,250 µl, transparent, PCR Performance**

**Tested, Low retention, 480 piece(s)/StackPack**



## Product information

Product number	70.3060.110
----------------	-------------

## Product descriptions

Name	Pipette tip, 1,250 µl, transparent, PCR Performance Tested, Low retention, 480 piece(s)/StackPack
Product number	70.3060.110
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	9 mm
Length of product	90 mm
Max. volume	1250 µl

Filter	no
Graduation	yes

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Piece(s) / inner box	1920
Piece(s) / case	1920
Piece(s) / pallet	184320
Height of case	225 mm
Width of case	267 mm
Depth of case	290 mm
Volume of case	0.0174 cbm
Weight of case	2.9702 kg
EAN inner box	4038917432238
EAN case	4038917432221

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT