




# Product specification

Pipette tip, 1,250 µl, transparent, PCR Performance

Tested, 500 piece(s)/bag



## Product information

Product number	70.3060
----------------	---------

## Product descriptions

Name	Pipette tip, 1,250 µl, transparent, PCR Performance Tested, 500 piece(s)/bag
Product number	70.3060
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 500 piece(s)/bag

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	9 mm
Length of product	90 mm
Max. volume	1250 µl

## Product characteristics

Filter	no
--------	----

Graduation	yes
Autoclavable (max. °C)	121 °C

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	120000
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	5.72 kg
EAN inner box	4038917434645
EAN case	4038917434638

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT