




Product specification

Filter tip, 1,000 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/box



Product information

Product number	70.3050.275
----------------	-------------

Product descriptions

Name	Filter tip, 1,000 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/box
Product number	70.3050.275
Product description	Filter tip, working volume: 1,000 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/box

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	9 mm
Length of product	72 mm
Max. volume	1000 µl

Product characteristics

Filter	yes
Graduation	yes

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	1920
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	61440
Height of inner box	148 mm
Width of inner box	111 mm
Depth of inner box	530 mm
Height of case	245 mm
Width of case	325 mm
Depth of case	547 mm
Volume of case	0.0436 cbm
Weight of case	5.1302 kg
EAN inner box	4038917424660
EAN case	4038917424653

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT