


# Product specification

**Pipette tip, 1000 µl, transparent, Biosphere® plus, 96 piece(s)/box**



## Product information

Order number	70.3050.205
--------------	-------------

## Product descriptions

Name	Pipette tip, 1000 µl, transparent, Biosphere® plus, 96 piece(s)/Box
Product number	70.3050.205
Product description	Pipette tip, working volume: 1000 µl, transparent, filling level rings, Biosphere® plus, suitable for Eppendorf, Gilson, Finnpiquette and Brand and products of identical design, 96 piece(s)/box

## Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C
Suitable for	Eppendorf, Gilson, Finnpiquette and Brand and products of identical design

## Size

Diameter	9 mm
Length of product	72 mm
Max. volume	1000 µl

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

## Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	61440
Height of case	245 mm
Width of case	325 mm
Depth of case	547 mm
Volume of case	0.0436 cbm
Weight of case	4.9901 kg
EAN inner box	4038917424752
EAN case	4038917424745

**This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.**

**This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.**

This document was prepared by EDP support and is valid without signature.

