




## Product specification

**Pipette tip, 1,000 µl, transparent, PCR Performance**

**Tested, Low retention, 480 piece(s)/StackPack**



## Product information

|                |             |
|----------------|-------------|
| Product number | 70.3050.110 |
|----------------|-------------|

## Product descriptions

|      |   |
|------|---|
| Name | Pipette tip, 1,000 µl, transparent, PCR Performance Tested, Low retention, 480 piece(s)/StackPack |
|------|---|

|                |             |
|----------------|-------------|
| Product number | 70.3050.110 |
|----------------|-------------|

|                     |  |
|---------------------|--|
| Product description | Pipette tip, working volume: 1,000 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack |
|---------------------|--|

## Material & colours

|                  |                    |
|------------------|--------------------|
| Product material | Polypropylene (PP) |
|------------------|--------------------|

|                  |             |
|------------------|-------------|
| Color of product | transparent |
|------------------|-------------|

## Size

|          |      |
|----------|------|
| Diameter | 9 mm |
|----------|------|

|                   |       |
|-------------------|-------|
| Length of product | 72 mm |
|-------------------|-------|

|             |         |
|-------------|---------|
| Max. volume | 1000 µl |
|-------------|---------|

## Product characteristics

|            |     |
|------------|-----|
| Filter     | no  |
| Graduation | yes |

## Purity & certification

|                  |                        |
|------------------|------------------------|
| Product category | diagnostic device, CE  |
| CE certified     | CE                     |
| Purity standard  | PCR Performance Tested |
| Batched          | yes                    |

## Packaging

|                        |               |
|------------------------|---------------|
| Minimum order quantity | 1920          |
| Piece(s) / inner box   | 1920          |
| Piece(s) / case        | 1920          |
| Piece(s) / pallet      | 149760        |
| Height of case         | 225 mm        |
| Width of case          | 295 mm        |
| Depth of case          | 303 mm        |
| Volume of case         | 0.0201 cbm    |
| Weight of case         | 2.7898 kg     |
| EAN inner box          | 4038917424813 |
| EAN case               | 4038917424806 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT