




Product specification

Pipette tip, 1,000 µl, transparent, PCR Performance

Tested, Low retention, 480 piece(s)/StackPack



Product information

Product number	70.3050.110
----------------	-------------

Product descriptions

Name	Pipette tip, 1,000 µl, transparent, PCR Performance Tested, Low retention, 480 piece(s)/StackPack
Product number	70.3050.110
Product description	Pipette tip, working volume: 1,000 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	9 mm
Length of product	72 mm
Max. volume	1000 µl

Filter	no
Graduation	yes

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	1920
Piece(s) / case	1920
Piece(s) / pallet	149760
Height of case	225 mm
Width of case	290 mm
Depth of case	303 mm
Volume of case	0.0198 cbm
Weight of case	2.7898 kg
EAN inner box	4038917424813
EAN case	4038917424806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT