




Product specification

Pipette tip, 1,000 µl, black, PCR Performance Tested, 500 piece(s)/bag



Product information

Product number	70.3050.030
----------------	-------------

Product descriptions

Name	Pipette tip, 1,000 µl, black, PCR Performance Tested, 500 piece(s)/bag
Product number	70.3050.030
Product description	Pipette tip, working volume: 1,000 µl, black, PCR Performance Tested, 500 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	black

Size

Diameter	9 mm
Length of product	72 mm
Max. volume	1000 µl

Product characteristics

Filter	no
Graduation	no
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	5000
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	180000
Height of case	300 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0475 cbm
Weight of case	5.24 kg
EAN inner box	4038917434614
EAN case	4038917434607

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT