




Product specification

Pipette tip, 1,000 µl, transparent, PCR Performance

Tested, 500 piece(s)/bag



Product information

Product number	70.3050
----------------	---------

Product descriptions

Name	Pipette tip, 1,000 µl, transparent, PCR Performance Tested, 500 piece(s)/bag
Product number	70.3050
Product description	Pipette tip, working volume: 1,000 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 500 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	9 mm
Length of product	72 mm
Max. volume	1000 µl

Product characteristics

Filter	no
--------	----

Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	180000
Height of case	310 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0491 cbm
Weight of case	4.76 kg
EAN inner box	4038917434553
EAN case	4038917434546

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT