




Product specification

Pipette tip, 300 µl, transparent, Biosphere® plus, 96 piece(s)/SingleRefill



Product information

Product number	70.3040.305
----------------	-------------

Product descriptions

Name	Pipette tip, 300 µl, transparent, Biosphere® plus, 96 piece(s)/SingleRefill
Product number	70.3040.305
Product description	Pipette tip, working volume: 300 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette, Biohit and Brand and products of identical design, 96 piece(s)/SingleRefill

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	7 mm
Length of product	60 mm
Max. volume	300 µl

Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	1920
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	103680
Height of inner box	148 mm
Width of inner box	105 mm
Depth of inner box	333 mm
Height of case	172 mm
Width of case	351 mm
Depth of case	437 mm
Volume of case	0.0264 cbm
Weight of case	2.7398 kg
EAN inner box	4038917433136
EAN case	4038917433129

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



SARSTEDT

This document was prepared by EDP support and is valid without signature.