




Product specification

Pipette tip, 300 µl, transparent, Biosphere® plus, 96 piece(s)/box



Product information

Product number	70.3040.205
----------------	-------------

Product descriptions

Name	Pipette tip, 300 µl, transparent, Biosphere® plus, 96 piece(s)/Box
Product number	70.3040.205
Product description	Pipette tip, working volume: 300 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	7 mm
Length of product	60 mm
Max. volume	300 µl

Product characteristics

Filter	no
--------	----

Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	61440
Height of inner box	148 mm
Width of inner box	111 mm
Depth of inner box	530 mm
Height of case	245 mm
Width of case	325 mm
Depth of case	547 mm
Volume of case	0.0436 cbm
Weight of case	4.0608 kg
EAN inner box	4038917433228
EAN case	4038917433211

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT