




Product specification

**Pipette tip, 300 µl, transparent, PCR Performance Tested,
Low retention, 480 piece(s)/StackPack**



Product information

Product number	70.3040.110
----------------	-------------

Product descriptions

Name	Pipette tip, 300 µl, transparent, PCR Performance Tested, Low retention, 480 piece(s)/StackPack
Product number	70.3040.110
Product description	Pipette tip, working volume: 300 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	7 mm
Length of product	60 mm
Max. volume	300 µl

Product characteristics

Filter	no
Graduation	yes

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	1920
Piece(s) / inner box	1920
Piece(s) / case	1920
Piece(s) / pallet	245760
Height of case	225 mm
Width of case	186 mm
Depth of case	290 mm
Volume of case	0.0121 cbm
Weight of case	1.6896 kg
EAN inner box	4038917433402
EAN case	4038917433396

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT