




Product specification

**Pipette tip, 300 µl, transparent, PCR Performance Tested,
1,000 piece(s)/bag**



Product information

Product number	70.3040
----------------	---------

Product descriptions

Name	Pipette tip, 300 µl, transparent, PCR Performance Tested, 1,000 piece(s)/Bag
Product number	70.3040
Product description	Pipette tip, working volume: 300 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 1,000 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	7 mm
Length of product	60 mm
Max. volume	300 µl

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	420000
Height of case	250 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0396 cbm
Weight of case	4.18 kg
EAN inner box	4038917434522
EAN case	4038917434515

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT