




Product specification

**Pipette tip, 200 µl, transparent, PCR Performance Tested,
1,000 piece(s)/bag**



Product information

Product number	70.3030
----------------	---------

Product descriptions

Name	Pipette tip, 200 µl, transparent, PCR Performance Tested, 1,000 piece(s)/bag
Product number	70.3030
Product description	Pipette tip, working volume: 200 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette, Biohit and Brand and products of identical design, 1,000 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	51 mm
Max. volume	200 µl

Product characteristics

Filter	no
--------	----

Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	480000
Height of case	205 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0301 cbm
Weight of case	3.83 kg
EAN inner box	4038917525947
EAN case	4038917525930

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT