



Product specification

Mailing container, transparent, construction: round, with absorbent liner, length: 114 mm, Ø opening



Product information

Product number	78.897
----------------	--------

Product descriptions

Name	Mailing container, transparent, construction: round, with absorbent liner, length: 114 mm, Ø opening: 44 mm, without cap
Product number	78.897
Product description	Mailing container, round, transparent, construction: round, with absorbent liner, length: 114 mm, Ø opening: 44 mm, without cap, cap: natural, 50 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Color of cap	natural

Size

Length of product	114 mm
Diameter of opening	44 mm

Product characteristics

Form of opening	round
Suction insert	yes
Example of use	Secondary packaging according to P 650
Liquid-tight	yes

Purity & certification

Satisfies the requirement	P 650, ADR
Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	250
Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	5000
Height of case	350 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0833 cbm
Weight of case	6.62 kg
EAN inner box	4038917084536
EAN case	4038917003346

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT