



Product specification

Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP, with print



Product information

Product number	62.547.274
----------------	------------

Product descriptions

Name	Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP, with print
------	---

Product number	62.547.274
----------------	------------

Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: natural, cap assembled, with print, label/print: blue/white, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/bag
---------------------	--

Material & colours

Color of product	transparent
------------------	-------------

Closure material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

Color of cap	natural
--------------	---------

Size

Diameter	28 mm
----------	-------

Length of product	114 mm
-------------------	--------

Base shape	Conical base
Label/Print	With print
Color of print/label	blue/white
Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	25
Piece(s) / case	300
Piece(s) / pallet	9000
Height of case	258 mm
Width of case	328 mm
Depth of case	580 mm
Volume of case	0.0491 cbm
Weight of case	4.61 kg
EAN inner box	4038917242561
EAN case	4038917242554

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT