



## Product specification

**Screw cap tube, 7 ml, (LxØ): 82 x 13 mm, PP**



### Product information

|                |            |
|----------------|------------|
| Product number | 60.550.050 |
|----------------|------------|

### Product descriptions

|                     |   |
|---------------------|---|
| Name                | Screw cap tube, 7 ml, (LxØ): 82 x 13 mm, PP   |
| Product number      | 60.550.050  |
| Product description | Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, flat base, transparent, screw cap, natural, cap enclosed, 1,000 piece(s)/bag |

### Material & colours

|                  |                                   |
|------------------|-----------------------------------|
| Color of product | transparent                       |
| Closure material | High Density Polyethylene (HD-PE) |
| Color of cap     | natural                           |

### Size

|                   |       |
|-------------------|-------|
| Diameter          | 13 mm |
| Length of product | 82 mm |

### Product characteristics

|                          |           |
|--------------------------|-----------|
| Base shape               | Flat base |
| Graduation               | no        |
| Closure type             | Screw cap |
| Centrifugation max (RCF) | 4000 x g  |

## Purity & certification

|                           |                         |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA                    |
| Product category          | no medical device   IVD |
| Batched                   | yes                     |

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 1000          |
| Piece(s) / case      | 1000          |
| Piece(s) / pallet    | 48000         |
| Height of case       | 215 mm        |
| Width of case        | 378 mm        |
| Depth of case        | 388 mm        |
| Volume of case       | 0.0315 cbm    |
| Weight of case       | 4.1 kg        |
| EAN inner box        | 4038917112536 |
| EAN case             | 4038917013734 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT