



## Product specification

**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP,  
without cap, 1,000 piece(s)/bag**



## Product information

Product number	60.504.021
----------------	------------

## Product descriptions

Name	Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag
Product number	60.504.021
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, without cap, 1,000 piece(s)/bag, 2,000 piece(s)/carton

## Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)

## Size

Diameter	13 mm
Length of product	75 mm

## Product characteristics

Base shape	Round base
------------	------------

Graduation	no
Closure type	Screw cap

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	72000
Volume of case	0.0492 cbm
Weight of case	5.74 kg
EAN inner box	4038917348201
EAN case	4038917348195

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**