



Product specification

Sample tube, EDTA K3E, 2 ml, cap red, (LxØ): 55 x 12 mm, with print



Product information

Product number	38.346
----------------	--------

Product descriptions

Name	Sample tube, EDTA K3E, 2 ml, cap red, (LxØ): 55 x 12 mm, with print
Product number	38.346
Product description	Sample tube, nominal volume: 2 ml, cap red, preparation: K3 EDTA, push cap, (LxØ): 55 x 12 mm, round base, with print, 100 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	red

Size

Diameter	12 mm
Length excluding cap	55 mm
Length of product	55 mm
Sample volume	2 ml

Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	red
Graduation	yes
Preparation concentration	1.6 mg/ml blood
Cap	Skirted push cap
Closure type	Push cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	2000
Piece(s) / pallet	72000
Height of case	265 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.042 cbm
Weight of case	5.81 kg
EAN inner box	4038917081283
EAN case	4038917000420

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

