



Product specification

Silicon adapter, for Automatic-Sarpette®



Product information

Product number	90.189.222
----------------	------------

Product descriptions

Name	Silicon adapter, for Automatic-Sarpette®
Product number	90.189.222
Product description	Silicon adapter, for Automatic-Sarpette®, translucent, silicone, 1 piece(s)/bag

Material & colours

Color of product	translucent
------------------	-------------

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / case	100
Piece(s) / pallet	32400
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm

Weight of case	1.65 kg
EAN inner box	4038917419086
EAN case	4038917419079

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT