



Product specification

Replacement filter, for Automatic-Sarpette®, pore size: 0.2 µm



Product information

Product number	90.189.221
----------------	------------

Product descriptions

Name	Replacement filter, for Automatic-Sarpette®, Pore size: 0.2 µm
Product number	90.189.221
Product description	Replacement filter, for Automatic-Sarpette®, pore size: 0.2 µm, 5 piece(s)/bag

Product characteristics

Pore size	0.2 µm
-----------	--------

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / case	500
Piece(s) / pallet	54000
Height of case	195 mm
Width of case	223 mm

Depth of case	308 mm
Volume of case	0.0134 cbm
Weight of case	1.205 kg
EAN inner box	4038917419604
EAN case	4038917419598

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT