

Product specification

Microvette® CB 300 Lithium heparin LH, 300 µl, cap orange, push cap, flat base



Product information

Product number	16.434
----------------	--------

Product descriptions

Name	Microvette® CB 300 Lithium heparin LH, 300 µl, cap orange, push cap, flat base
------	--

Product number	16.434
----------------	--------

Product description	Microvette® CB 300 Lithium heparin LH, Veterinary medicine, preparation: Lithium heparin, nominal volume: 300 µl, (LxØ) with cap: 46.5 x 10.8 mm, cap orange, push cap, flat base, inner receptacle: trumpet-shaped, carrier tube: brown, 100 piece(s)/bag
---------------------	--

Material & colours

Color of product	transparent
------------------	-------------

Closure material	Low Density Polyethylene (LD-PE)
------------------	----------------------------------

Color of cap	orange
--------------	--------

Size

Diameter	10.8 mm
----------	---------

Length excluding cap	43.7 mm
----------------------	---------

Length including cap	46.5 mm
----------------------	---------

Sample volume 300 µl

Product characteristics

Base shape	Flat base
Shape of inner tube	trumpet-shaped
Cap	Soft stopper with retaining loop
Closure type	Push cap
Type of collection	capillary
Centrifugation recommended	2000 U/min
Centrifugation time	5 Min

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	80000
Height of case	180 mm
Width of case	288 mm
Depth of case	378 mm
Volume of case	0.0196 cbm
Weight of case	2.35 kg
EAN inner box	4038917268264
EAN case	4038917282697

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

