



Product specification

M-Sarpette®, Bulb Pipettor



Product information

Product number	90.1427
----------------	---------

Product descriptions

Name	M-Sarpette®, Bulb Pipettor
Product number	90.1427
Product description	M-Sarpette®, Bulb Pipettor, 1 piece(s)/case

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / case	40
Piece(s) / pallet	1440
Height of case	300 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0475 cbm
Weight of case	6.69 kg
EAN inner box	4038917086592
EAN case	4038917005135

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT