




## Product specification

### Stripping / Sealing Tongs



#### Product information

Product number	97.8800.200
----------------	-------------

#### Product descriptions

Name	Stripping / Sealing Tongs
Product number	97.8800.200
Product description	Stripping / Sealing Tongs, for stripping, segmenting and sealing blood collection tubes

#### Material & colours

Color of product	red
------------------	-----

#### Size

Width of product	5 cm
Length of product	16 cm

#### Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

#### Packaging

Piece(s) / pallet	432
Height of case	100 mm

Width of case	188 mm
Depth of case	188 mm
Volume of case	0.0035 cbm
Weight of case	0.264 kg
EAN inner box	4038917975988
EAN case	4038917975995

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT