




Product specification

Sealing clips



Product information

Product number	88.7708.095
----------------	-------------

Product descriptions

Name	Sealing clips
Product number	88.7708.095
Product description	Sealing clips, for sealing blood collection tubes, material: aluminium

Material & colours

Color of product	grey
------------------	------

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / inner box	1000
Piece(s) / case	99000
Piece(s) / pallet	6336000
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm

Weight of case	11.88 kg
EAN inner box	4038917140454
EAN case	4038917550222

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT