




Product specification

Tube stripper AST-2F with foot switch



Product information

Product number	87.7709.035
----------------	-------------

Product descriptions

Name	Tube stripper AST-2F with Foot switch
Product number	87.7709.035
Product description	Tube stripper AST-2F, with auto-stop function, foot switch, table holder

Size

Width of product	140 mm
Length of product	250 mm
Height of product	71 mm

Product characteristics

Power consumption	20 W
Supply voltage	240 V

Purity & certification

Product category	no medical device IVD
CE certified	CE
Batched	no

Packaging

Piece(s) / pallet	84
-------------------	----

Height of case	125 mm
Width of case	358 mm
Depth of case	388 mm
Volume of case	0.0174 cbm
Weight of case	2.76 kg
EAN inner box	4038917258890
EAN case	4038917258883

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT